

## TriM Scope

### Beam Multiplexer for Multifocal Multiphoton Microscopy

#### description

LaVision BioTec's TriMScope is based on a patented beamsplitter that splits up an incoming laser beam into up to 64 beams which are scanned simultaneously in the object plane. This results in either 64 times brighter images or 64 times higher image acquisition rates compared to standard single beam multiphoton scanning microscopes. The foci in the object plane are aligned in a single line and the number of foci can be easily switched from 64 to 32, 16, 8, 4 and to a single beam. With optional resonant scanner image rates up to 4000 Hz are possible.

The TriMScope beam divider is a compact and easy-to-use device that utilizes exclusively flat optics for dividing the incoming beam with high light efficiency, avoiding aberrations and producing equally intense foci in the sample. As the beams are inherently shifted with respect to each other by several picoseconds, there is no cross-talk. Spatial and axial resolutions are diffraction limited. In addition, single line excitation can easily be coupled into the input slit of an imaging spectrograph yielding real-time spectral sectioning in the x-y and x-z planes.

#### microscopes

invert and upright microscopes from Leica, Nikon, Olympus and Zeiss

#### laser

Ti:Sapphire lasers from Spectra-Physics, Coherent, etc.

#### frame rate

The frame rate depends on the CCD camera.

**frame rate:** with CCD-camera Andor iXon DV885  
31 Hz @ 1004 x 1002 pixel

#### scan rate

**max. scan rate:** 1.2 kHz with non resonant scanner

#### scan area/ FOV

The FOV depends on the CCD chip size and the magnification of the objective. The following values correspond to the Andor iXon DV 885 CCD camera (chip size 8.0 x 7.0 mm<sup>2</sup>).

**100 x Objective:** 80 x 80 μm<sup>2</sup>  
**63 x Objective:** 127 x 127 μm<sup>2</sup>  
**20 x Objective:** 400 x 400 μm<sup>2</sup>

#### transmission

> 70% within the TriM Scope  
@ 710 – 980 nm

#### z-drive

objective driven via a precision stepper motor

**speed:** 10 Hz @ 0.9 μm step size  
**min. step size:** 25 nm

#### time multiplexing

**time delay  
between adjacent  
beams** 680 fs

<b>inter-foci separation</b>	<p>100 x Objective: 300 nm  63 x Objective: 476 nm  20 x Objective: 1500 nm</p> <p>The focal spot size is diffraction limited and depends on the objective.</p>
<b>foci size</b>	<p>63 x Objective: <math>\approx 300</math> nm  @ <math>\lambda_{\text{laser}} = 800</math> nm  NA<sub>objective</sub>: 1.4</p>
<b>polarization of adjacent foci</b>	<p>The polarization of foci alternates between S and P.</p> <p>Arrangement: s-p-s-p-s-p...</p>
<b>variable number of foci</b>	<p>The foci are arranged in a line. The number of foci can be set via software. The length of the line is proportional to the number of foci.</p> <p>number of foci: 1, 2, 4, 8, 16, 32, 64  point/line scan toggle switch</p>
<b>chirp compensation</b>	<p>optical dispersion causes broadening of the laser pulse leading to a linear decrease in 2-photon induced fluorescence signal. TriMScope has a chirp compensation arrangement based on prisms.</p> <p>chirp compensation: laser pulse width down to 120 fs at the sample plane based on the original pulse width of 100 fs.</p>
<b>intensity control</b>	<p>via motorized crossed polarizer arrangement at the input</p>
<b>power requirements</b>	<p>110V/220V ; 5A</p>
<b>dimensions</b>	<p><math>60 \times 42 \times 18</math> cm<sup>3</sup></p>
<b>weight</b>	<p>42 kg</p>
<b>upgrade spectral unmixing</b>	<p>The line of 64 foci lends itself for convenient coupling onto the input slit of an imaging spectrograph allowing spectral-resolved imaging.</p>
<b>upgrade FRAP</b>	<p>An additional visible laser can be coupled in the optical path which can be scanned via the XY scanner and used for FRAP studies or uncaging.</p>
<b>upgrade 3D FLIM</b>	<p>TriM Scope in conjunction with LaVision's ultrahigh rep. rate, picosecond gated PicoStar HR (200ps gate width @ 110 MHz rep. rate) camera allows real-time 3D FLIM.</p>
<b>customized systems</b>	<p>We specialize in offering customized turn-key systems. Please feel welcome to contact us to discuss your special needs and requirements.</p>